



EASY, FAST, RELIABLE

ELDONCARD® 2511 - ABO & RhD forward blood typing



Recommended for:

- Medical Market – Blood Banks, Hospitals, Clinics etc.
- Emergency & Rescue – NGO's
- Military, Police and Civil Service
- Donor Recruitment
- Rhesus tests for Pregnant Women
- Educational Market
- Health / Home Market
- Bed side test

The card is used for the detection of the blood type of one person
Control field to detect false possible false positives.

The procedure is described in detail in the "Instruction for Use".

Instructional video: www.eldoncard.com (English/Espanol/French + prof. use)

Available as:



- **ELDON BAG**
Contains:
 - 25 ELDONCARD 2511
25 Blood Group Determinations
- **ELDON BOX 2511-25**
Contains:*
 - 100 ELDONCARDS 2511
100 Blood Group Determinations
- **ELDON BIG BAG**
Contains:*
 - 50 or 100 ELDONCARDS 2511
50/100 Blood Group Determinations

- **Costume made Kits for special use NGO/Doctors/ Military**

**Each EldonCard is in an aluminum envelope.*

Labelling: CE marked & FDA cleared + a number of countries

Storage: 5-37 Celsius / 41-99 Fahrenheit
NO COOL CHAIN is needed

Expiration: 24 month from first day of production

Accessories:

In general, accessories are available in numbers matching the number of Blood Group Determinations.



- Safety Lancet



- Steel Lancet



- EldonStick



- Alcohol Swabs



- EldonFoil



- Pipette

General: To perform one Blood Group Determination 2511 you need 1 lancet, 4 EldonSticks, 1 EldonFoil, 1 pipette and 1 alcohol swab.