


ELDONCARD® 2521 - for ABO & RhD forward blood typing



Recommended for:

- ☑ Military, Police and Civil Service
- ☑ Donor Recruitment
- ☑ Rhesus tests for Pregnant Women
- ☑ Educational Market
- ☑ Health / Home Market
- ☑ Bed side test

ELDONCARD® 2521 carries 2 x 4 circles with a diameter of 20 mm.

The card is used for the detection of the blood type of one/two person/s Control field to detect possible false positives.

The procedure is described in detail in the "Instruction for Use".

Instructional video: www.eldoncard.com (English/Espanol/French + prof. use)

Available as:



- ELDON BAG 2521 Contains:
 - 25 ELDONCARD 2521 25 Blood Group Determinations
- ELDON BOX 2521 Contains:*
 - 100 ELDONCARDS 2521
 100 Blood Group Determinations





*Each EldonCard is in an aluminum envelope.

Labelling: CE marked & FDA cleared and in other countries **Storage:** 5-37 Celsius / 41-99 Fahrenheit – **NO COOL CHAIN Expiration:** 24 month from first day of production

Accessories:

In general, accessories are available in numbers matching the number of Blood Group Determinations.



Safety Lancet



Steel Lancet



EldonStick



Alcohol Swabs



EldonFoil



Pippette



General: To perform one Blood Group Determination 2521 you need 1+1 lancet, 4x2 EldonSticks, 1 EldonFoil, 1 pipette and 1+1 alcohol swab.