



EASY, FAST, RELIABLE

ELDONCARD® RhD: Forward blood typing



Recommended for:

- Medical Market – Blood Banks, Hospitals, Clinics etc.
- Emergency & Rescue – NGO's
- Military, Police and Civil Service
- Donor Recruitment
- Rhesus tests for Pregnant Women
- Educational Market
- Health / Home Market
- Bed side test

ELDONCARD® RhD. Each pair of fields contain Anti-D and Control reagents, respectively. The card has two perforation lines so that it easily can be divided into three individual cards and is intended for blood typing of three individuals. The intended use is RhD determination of 3 individuals per card. Control field to detect possible false positives.

Available as:



- **ELDON BAG Contains:***
- 25 ELDONCARD RhD
25 Blood Group Determinations
- BIG BAG ELDONCARD RhD
100 Blood Group Determinations

- **Costume made Kits for special use NGO/Doctors/ Military**

**Each EldonCard is in an aluminum envelope.*

Labelling: CE marked & FDA cleared + a large number of other countries

Storage: 5-37 Celsius / 41-99 Fahrenheit
NO COOL CHAIN

Expiration: 24 month from first day of production

Accessories:

In general, accessories are available in numbers matching the number of Blood Group Determinations.



- Safety Lancet



- Steel Lancet



- EldonStick



- Alcohol Swabs



- EldonFoil



- Pipette

- The procedure is described in detail in the "Instruction for Use".
- Instructional video: www.eldoncard.com (English/Espanol/French + prof. use)